

CASE STUDY

CCI Test Method Validation for Pre-filled Syringes





Background

- Regulators are paying closer attention to the proper design of robust Container Closure Integrity (CCI) studies and the validation of CCI test methods.
- A method was developed for CCI testing of pre-filled syringes and then validated.
- The purpose of Method Validation is to demonstrate sufficient accuracy and reproducibility of a method, in this case to detect CCI defects at the plunger seal of a pre-filled syringe.
- The developed CCI test method uses Headspace Gas Ingress analysis to identify changes in the headspace content of the syringe that result from gas ingress through a leak.



The Method

The method applied uses an FMS Headspace Gas Analyzer and a CCI Test Vessel, offering a robust solution for CCI Testing. This method utilizes Headspace Gas Ingress analysis to measure changes in the headspace gas composition of the syringe that result from gas ingress through a leak.

The CCI Test Vessel creates an environment that accelerates gas ingress through an open defect. Samples are placed inside the vessel and exposed to a tracer gas (in this case carbon dioxide). In the case of a closure defect, CO₂ will leak into the syringe headspace and will be detected by performing headspace analysis.

METHOD VALIDATION

Three operators performed headspace analysis, each on 5 samples of every control group.

Before measurement, the samples were pre-conditioned in the CCI Test Vessel for 30 minutes at 1.3 atm of CO₂ overpressure and subsequent headspace analysis of all the samples was easily performed in less than 15 minutes.

Sample overview:

NUMBER OF SAMPLES	CONTROL GROUPS
15	Positive Control - small defect - empty (Figure 1)
15	Positive Control - small defect - formulation
15	Positive Control - small defect - empty
15	Negative Control



Figure 1
Positive control with a tungsten microwire creating a small defect.






Conclusions

The purpose of the Method Validation was to demonstrate that the method is sufficient for accurately and reproducibly detecting critical defects at the plunger seal. All the positive controls in this study measured more than half an atmosphere of CO₂ ingress into the syringe headspace. All negative controls measured no CO₂ ingress.

This headspace gas ingress CCI testing method was shown to be validated for the detection of a critical leak for a specific product configuration. The parameters of the method can be straightforwardly optimized to generate robust CCI data on other products if appropriate.



1 ML SYRINGES	CO ₂ INGRESS (TORR)	LEAK	CCIT
	AVERAGE	DETECTED	RESULT
Positive Control - empty	451	15/15	100% Pass
Positive Control - filled	573	15/15	100% Pass
Gross Defect	690	15/15	100% Pass
Negative Control	0	0/15	100% Pass

CONTACT US

- Lighthouse Instruments, LLC
 2050 Avon Court
 Charlottesville, VA 22902
 United States of America
Office • +1 (434) 293-3081
Email • info@lighthouseinstruments.com

- Lighthouse Instruments B.V.
 Wisselwerking 22
 1112 XP Diemen
 The Netherlands
Office • +31 (0)20 7051 050
Email • euinfo@lighthouseinstruments.com

www.lighthouseinstruments.com