


PRODUCT NOTE

FMS-WATER ACTIVITY Headspace Analyzer





The FMS-Water Activity Analyzer is a non-destructive gas analyzer for making water activity measurements of solid dosage drug product samples. This compact analyzer utilizes a patented laser absorption technique developed with funding from the Food and Drug Administration. Light from a near-infrared laser is tuned to match an internal absorption frequency of the water vapor molecule and passed through a glass sample container in the headspace above the product. The amount of laser light absorbed is proportional to the water vapor concentration in the headspace which in turn can be directly related to the water activity of the pharmaceutical solid. The non-destructive nature of the measurement method enables a broad range of moisture activity applications and offers significant advantages over traditional techniques.

NOMINAL SPECIFICATIONS	
Moisture Measurement Linearity Range	0 to 100% RH
Temperature Control Range	19 to 45°C
Measurement Time	5 seconds (default)
Sample container sizes	1 ml to 50 ml
PHYSICAL AND ELECTRICAL	
Dimensions (H×W×D)	292 × 305 × 405 mm (11.5" × 12.0" × 15.9")
Weight	19.5 kg (43 lbs)
Power Requirements	110–240 VAC, 50/60 Hz, 180W
Controller	PC

APPLICATIONS INCLUDE:

- Water activity determination of solid dosage and powder product
- Real-time moisture exchange studies
- Water vapor transmission rate studies
- Replacement of Karl Fischer Titration (KFT), Loss On Drying (LOD) for lyo moisture determination
- Rubber stopper drying studies



PRODUCT NOTE

FMS-Water Activity Headspace Analyzer



STANDARDS	ACTUAL (%RH)	MEAN (%RH)	STEV (%RH)
water	100	98.8	0.28
KN03	93.6	93.3	0.78
NaCl	75.3	76.8	1.56
NaBr	57.6	59.9	1.20
MgCl2	32.8	33.4	0.99
LiCl	11.3	11.6	0.07
LiBr	6.4	6.8	0.14
Desiccant	0.0	0.6	0.42



Table 1. Relative humidity measurements on known moisture standards

MOISTURE PERFORMANCE DATA

Six clear glass vials were prepared with saturated salt solutions to create moisture standards at approximately 94, 75, 58, 33, 11, and 6% RH. Water and desiccant samples were used for 100% and 0% RH standards. The relative humidity in each vial was measured at an equilibration temperature of 24 degrees Celsius using the FMS-Water Activity Analyzer. The analyzer measures the absorption of laser light by water vapor molecules in the container headspace. The amplitude of the absorption signal is proportional to the relative humidity. The mean measured moisture for each of the eight moisture standards is presented in Table 1.

Figure 1 is a plot of the measured versus known relative humidity in each vial. The non-destructive nature of the measurement enables a full statistical documentation of system performance. From an application point-of-view, the non-destructive measurement enables multiple measurements on the same sample over time, 100% analysis of a batch giving insight into process variability and optimization, and real-time monitoring of moisture dynamics.




Figure 1. Linearity for moisture measurements by the FMS-Water Activity Analyzer.





CONTACT US

- Lighthouse Instruments, LLC
2050 Avon Court
Charlottesville, Virginia 22902
United States of America
Office • +1 (434) 293-3081
Email • info@lighthouseinstruments.com

- Lighthouse Instruments, B.V.
Wisselwerking 22
1112 XP Diemen
The Netherlands
Office • +31 (0)20 7051 050
Email • euinfo@lighthouseinstruments.com

www.lighthouseinstruments.com

